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
INSTRUCTIONS FOR USE FOR VISUALLY READ SERVOTESTURINE TEST STRIPS

USAGE IN VITRO CE ⊗ IVD

SUMMARY AND EXPLANATION

Reagent strips for the rapid determination of Urobilinogen, Glucose, Bilirubin, Ketones, Specific Gravity, Blood, pH, Nitrite, Leucocytes and Ascorbic Acid. The combination of the different parameters are shown on the vial label/packaging.

INTENDED USE

SERVOTEST Reagent Strips are dip-and-read test strips for In Vitro Diagnostic Use only for testing above items in urine. Test result may provide information regarding the status of carbohydrate metabolism, kidney and liver function, acid-base balance, and urinary tract infection. It is measured by comparison of test paper attached to a plastic strip with the color chart blocks printed on the vial label. The strips may be read visually.

STORAGE AND SHELF LIFE

- Replace the bottle cap immediately and tightly after removing test strips, and keep the vial tightly closed between tests.
- Store in a cool, dry place at temperatures between 15 - 30 °C. Do not store the strips in a refrigerator or freezer.
- Store away from moisture and sunlight.
- When stored in the original container, the product is stable up to the expiry date printed on the bottom of the container
- Do not use after expiration date.

Discoloration or darkening of the test pads may indicate deterioration. If this is evident, or if test results are questionable or inconsistent with expected finding, confirm that the product is within its expiration date and is reacting properly using known negative and positive control materials.

LIMITATIONS OF PROCEDURE

As with all laboratory tests, definitive diagnostic or therapeutic decisions should not be based on any single result of method. Substances that cause abnormal urine color may affect the readability of test pads in urinalysis reagent strips. Urinary ascorbic acid concentrations as low as 50mg/dl can cause interference in specimens with low concentrations of glucose, blood and bilirubin.

ANALYTICAL PRECAUTIONS

- Reagent strips are for diagnostic use only and should not be used for the analysis of body fluids other than urine.
- As with all laboratory tests, definitive diagnostic or therapeutic decisions should not be based on any single result of method.
- The effects of drugs or other metabolites on the individual tests of the test strips are not known in all cases. It is therefore recommended that in case of doubt, the test should be repeated after withdrawal of the potential interfering agent such as medication or vitamin supplement, etc.
- Do not remove desiccant from bottle.
- Do not touch test areas of urine reagent strips.
- Do not open container until ready to use.
- The work area should be clean and free from detergents and other contaminants.
- Keep out of reach of children.
- Each test strip is for single use only.
- The correct reading time shown on the vial label is important for optimal results and readings outside this will invalidate the test.
- Colour changes that appear only along the edge of the test area should be ignored careful removal of excess urine should eliminate this phenomenon.

SPECIMEN COLLECTION AND PREPARATION

Collect urine in a clean, dry container that allows complete immersion of all the fields on the test strip. Do not add preservatives. Test the specimen as soon as possible, with the sample well mixed but not centrifuged. The use of fresh morning urine is recommended for optimal nitrite tests, as well as for the valid determination of bilirubin and urobilinogen, since these compounds are unstable when exposed to light. If immediate testing is not possible, the sample should be stored in the refrigerator, but not frozen, and then brought to room temperature before used in the test. Unpreserved urine at room temperature may undergo pH changes due to microbial proliferation, which may interfere with protein determination. If cleanly voided specimens are not collected from females, positive results for leukocytes may be found due to contamination from outside the urinary tract. Skin cleansers containing chlorhexidine may affect protein test results if specimen contamination occurs.

HEALTH AND SAFETY WARNINGS

- All patient samples should be treated as potentially infectious and the user should wear appropriate protective equipment when performing the test.
- The reagents which are impregnated into each pad together with average quantities are listed in each section describing the principles of each section describing the principles of each test.

WARNING AND PRECAUTIONS

- For in vitro diagnostic use only.
- For professional use only.


QUALITY CONTROL

For best results, performance of reagent strips should be confirmed by testing known negative and positive specimen or controls whenever a new bottle is first opened. Each laboratory should establish its own goals for adequate standards of performance. Each lab worker should ensure that it complies with government and local requirements.

VISUAL TEST PROCEDURE

The procedure must be followed exactly to achieve reliable results.

1. Remove a strip from the vial and replace the cap immediately; Inspect the strip. If reagent areas are discolored or darkened, do not use the strip.
2. Refer to illustrations (1, 2, 3) for the following steps.
 - (a) Dip the strip into the urine up to the test area for no more than two second.
 - (b) Draw the edge of the strip along the brim of the vessel to remove excess urine; At this time, don't make the test areas touched to the brim of the vessel.
 - (c) Turn the strip on its side and tap once on a piece of absorbent material to remove any remaining urine; Excessive urine on the strip may cause the interaction of chemicals between adjacent reagent pads, so that an incorrect result may occur.
 - (d) Compare the colours of the reagent pads within 60 seconds (Leucocytes after 120 seconds) with the colour chart on the vial under good light. While comparing, keep the strip horizontally to prevent possible mixing of chemicals when excessive urine is present.

REAGENT AREA INFORMATION

1. UROBILINOGEN

Chemical Principle: The test is based on the Ehrlich's reaction. Colour changes from light orange-pink to dark pink.

Ingredients: 4-Methoxybenzenediazonium 2.9mg

Expected Values: The normal urobilinogen range is 0.1 to 1.0 Ehrlich unit /dl. If

results exceed the concentration of 2.0 mg/dl, the patient and the urine specimen should be evaluated further.

Limitations: The absence of urobilinogen in the specimen cannot be determined. The test area will react with interfering substances known to react with Ehrlich's reagent, such as p-aminosalicylic acid. Drugs containing azo gantrisin may give a masking golden color. The test is not reliable method for the detection of porphobilinogen.

2. GLUCOSE

Chemical Principle: Glucose oxidase catalyzes the oxidation of glucose to form hydrogen peroxide. The hydrogen peroxide thus formed then oxidizes a chromogen on the reaction pad by the action of peroxidase.

Ingredients: Glucose oxidase 430 U, Peroxidase 200 U, Potassium Iodide 12 mg

Expected Values: The kidney normally excretes small amounts of glucose. Concentrations of 100 mg/dl may be considered as abnormal if found consistently.

Limitations: High SG (>1.020) with high pH urine and ascorbic acid (more than 40 mg/dl) may cause false negative result at the low level of glucose. Ketones reduce the sensitivity of the test. Moderately high ketone level (>40 mg/dl) may cause a false negative for specimen containing small amount of glucose (100 mg/dl). Reactivity may be influenced by urine SG and temperature.

3. BILIRUBIN

Chemical Principle: Azo-coupling reaction of bilirubin with a diazonium salt in an acid medium to form an azodye. Color changes from light tan to beige or light pink.

Ingredients: Sodium nitrite 0.733 mg, 2,4-dichlorobenzene diazonium 2.3 mg, Sulfosalicylic acid 25 mg.

Expected Values: Normally no bilirubin is detectable. Even trace amounts of bilirubin are sufficiently abnormal to require further investigation.

Limitations: Metabolites of drugs, such as pyridium and serenium, which give a color at low pH, may cause false positives. Indican indoxyl sulfate can produce a yellow-orange to red colour response, which may interfere with the interpretation of negative or positive bilirubin readings. Ascorbic acid concentrations (>25 mg/dl) may cause false negative results.

4. KETONES

Chemical Principle: Legal's test-nitroprusside reaction. Acetoacetic acid in an alkaline medium reacts with nitroferric anion to produce a colour change from beige to purple.

Ingredients: Sodium nitroprusside 23.0 mg

Expected Values: Ketone bodies should not be detected in normal urine specimens with this reagent.

Limitations: Positive results (trace or less) may occur with highly pigmented urine specimens or those containing large amounts of levodopa metabolites. Some high SG and low pH urine may give false positive result. Phenosulfonphthalein may cause false positive result.

5. PH

Chemical Principle: Double indicator system. Indicator's methyl red and bromothymol blue are used to give distinct colour changes from orange to green to blue. (pH 5.0 to 9.0)

Ingredients: Methyl red 0.05 mg, Bromothymol blue 0.5 mg

Expected Values: Urine values generally range from pH 5 to 9.

Limitations: If the excessive urine is remain on the strip because of improper test procedure, it is possible that the acidic buffer in protein portion comes out and affect the pH portion, then pH result may be decreased than the actual. This phenomenon is called "run-over effect".

6. BLOOD

Chemical Principle: The test is based on the Pseudo-peroxidase activity of the haem moiety of hemoglobin and myoglobin. The chromogen is oxidized by a hydroperoxide in the presence of haem and changes colour from yellow to blue.

Ingredients: Cumene Hydroperoxide 12 mg, o-Tolidine 35 mg

Expected Values: Normally, no hemoglobin is detectable in urine (0.010 mg/dl; 3 RBC/ μ l). When hemoglobin appears in urine it indicates kidney disease or a urinary tract disorder. Blood may often be found in the urine of menstruating females.

Limitations: Elevated specific gravity or protein in urine may reduce the reactivity of the blood test portion. Microbial peroxidase associated with urinary tract infection may cause false positive results. Ascorbic acid concentrations (>40 mg/dl) may cause false negatives at the low level of blood.

7. SPECIFIC GRAVITY (SG)

Chemical Principle: Ionic solutes present in the urine cause protons to be released from a polyelectrolyte. As the protons are released the pH decreases and produces a colour change of bromothymol blue from blue-green to yellow-green.

Ingredients: Bromothymol blue 0.5 mg, Poly vinyl ether-ALT-maleic acid anhydrous 140.5 mg

Expected Values: The normal SG of urine ranges from 1.001 to 1.035

Limitations: High-buffered alkaline urine may cause diminished result, whereas high-buffered acidic urine may cause slightly elevated result.

8. PROTEIN

Chemical Principle: Protein "error of indicators." When pH is held constant by a buffer, indicator dyes release H⁺ ions because of the protein present and change colour from yellow to blue-green.

Ingredients: Tetrabromophenol blue 0.34 mg

Expected Values: Normal urine specimens ordinarily contain some protein (<20 mg/dL) therefore only persistent elevated levels of urine protein indicate kidney or urinary tract disease. The persistent results of trace level or over indicate significance proteinuria and thus further clinical testing is needed to evaluate the significant of results.

Limitations: False positive results may be found in strongly basic urine (pH 9). The interpretation of results is also difficult in turbid urine specimens.

9. NITRITE

Chemical Principle: The test is based on the diazotization reaction of nitrite with an aromatic amine to produce a diazonium salt. It is followed by an azo-coupling reaction of this diazonium salt with an aromatic compound on the reaction pad. The azo-dye produced causes a colour change from white to pink.

Ingredients: P-arsanilic acid 4.5 mg

Expected Values: Normally no nitrite is detectable in urine.

Limitations: Ascorbic acid (>25 mg/dL) may cause false negative result with low level of nitrite containing (<0.03 mg) urine. The negative result does not always mean that the patient is free from bacteriuria. Pink spots or pink edges should not be interpreted as a positive result. Negative result may occur when urinary tract infections are caused by organism which do not contain nitrate reductase; when urine has not been retained in the bladder long enough (four hours or more) for reduction of nitrate to nitrite occur; or when dietary nitrate is absent.

10. LEUKOCYTE

Chemical Principle: This test pad contains an indoxyl ester and diazonium salt. It is followed by an azo-coupling reaction of the aromatic amine formed by leukocytes esterase with a diazonium salt on the reaction pad. The azo dye produced causes a colour change from beige to violet.

Ingredients: Induced indole amino acid ester 1.3 mg

Expected Values: Normally no leukocytes are detectable in urine.

Limitations: The test result may not always be consistent with the leukocyte cell number by the microscopic examination. High concentration of glucose, high specific gravity, high level of albumin, high concentration of formaldehyde or presence of blood may cause decreased test results. High concentration of oxalic acid of trace of oxidizing agents may cause false negative results.

11. ASCORBIC ACID

Chemical Principle: The test field involves the decolorization of Tillmann's reagent. The presence of ascorbic acid causes the color of the test field to change from gray-blue to yellow.

Ingredients: 2,6-dichloro indophenol sodium salt 0.8 mg

Expected Values: The average daily intake ranges from 30-80 mg, with an output of 20-30 mg/day.

Limitations: No interferences are known.

PERFORMANCE CHARACTERISTICS

Performance characteristics are based on clinical and analytical studies and depend upon several factors: the variability of colour perception; the presence or absence of inhibitory and matrix factors typically found in urine; and the laboratory conditions in which the product is used (e.g., lighting, temperature, and humidity). Each colour block represents a range of values. Because of specimen and reading variability, specimens with analyte concentrations that fall between normal levels may give results at either level. Results will usually be within one level of the true concentration. The following list shows the generally detectable levels of the analytes in contrived urines; however, because of the inherent variability of clinical urines, lesser concentrations may be detected under certain conditions.

SERVOTEST TEST PAD AND SENSITIVITY (SPECIFICITY)

Urobilinogen:	2 IU/dL	(Urobilinogen)
Glucose:	100 mg/dL	(Glucose)
Bilirubin:	1 mg/dL	(Bilirubin)
Ketones:	5 mg/dL	(Acetoacetic acid)
Blood:	10 RBC/ μ l	(0.03 mg/dL hemoglobin, Intact RBC)
Protein:	10 - 15 mg/dL	(Albumin)
Nitrite:	0.05 mg/dL	(Nitrite ion)
Leukocytes:	20 - 25 WBC/ μ l	(Intact and lysed WBCs)
Ascorbic acid:	20 mg/dL	(Ascorbic acid)

Index of symbols

	Article number		Temperature limitation
	Observe operating instructions		Batch number
	In-vitro-diagnostic		Expiry date
	Manufacturer		Content sufficient for <n> tests
	Dangerous substances		Single use
	Protect from heat and sunlight		Attention
	Protect from moisture		
	Do not use, if package is damaged		
	CE marked according to IVD directive 98/79/EG		

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servoprax GmbH