



Supertronic

SUPERTRONIC IS A
DIAGNOSTIC DEVICE ACCORDING TO THE EAV
(ELECTROACUPUNCTURE ACCORDING TO DR. VOLL)



EAV - energy diagnosis

Electroacupuncture according to Voll is the oldest and most frequently used electrical diagnosis and therapy method in naturopathy with which acupuncture points can be precisely located and their energy state measured. It was developed according to the basic idea of Traditional Chinese Medicine: The precisely defined points on the surface of the skin, called acupuncture points, are connected to organs via energy pathways (meridians). There is therefore a reflexive relationship or interaction between the two.

The EAV device measures the conductivity of the tissue (organ resistance, in these special points). The measured value allows diagnostic conclusions to be drawn about the functional state of the corresponding organ (organ system).



Dr. REINHOLD VOLL

Dr. Reinhold Voll (1909 - 1989) was a Bavarian physician who was the first expert to combine the methods of conventional and unconventional medicine into a sophisticated whole. He used medical electrophysiological methods for acupuncture system research and defined over 800 points, with which the associated internal organs and tissues were connected.

The measurements allow us to determine if the organs, functional systems or tissues are irritated, inflamed, adjusting normally or degenerating. Energy diagnostics from the perspective of Voll's method examines the acupuncture point response to the direct current passage with a very low voltampere characteristic. This provides a representation of both the energy of the acupuncture pathway and the condition of the associated organs and tissues.



Häufig gestellte Fragen

When is a Supertronic examination most useful?

This method is most beneficial for prevention and in cases where conventional medicine does not identify the cause of a disease or does not know how to cure it.

Are conventional and unconventional medicine contradictory or complementary?

Conventional and unconventional medicine do not compete with each other, nor are they at odds. The fact that they work from different angles means that they actually expand the options for diagnosing and treating a disease. When used together, this is ideal for the person being studied, as the causes of the pathological condition can be identified more effectively and tailored treatment can be proposed.

How should you prepare for a Supertronic examination?

Patients should not be physically or mentally exhausted or suffer from an infectious disease. Do not smoke or drink for at least two hours before the examination and do not consume caffeine or energy drinks. Do not apply any cream to your hands or feet. Take only essential medications on the day of the examination. Do not take any other medicines or supplements. The counselor or therapist with whom you make your appointment will give you more details.

Are there any limiting factors or people?

Limiting factors for the examination are age (this depends on the size of the fingers of small children), pathological processes visible on the skin at the sites where measurements are to be taken (eczema, skin rashes or skin lesions or scars) and pregnancy (which could influence the interpretation of the measurements). Patients with pacemakers or other electronic devices or patients suffering from epilepsy should not be tested.



Supertronic APP

- » iOS - application for download in App Store in English language.
- » database of customers
- » individual settings of personal profile of measurement points
- » automatic reading and recording of measurements
- » synoptic analysis of results
- » comparison with previous measurements
- » recording of recommended products
- » user-friendly interface



Compatibility:

iPad 3 and new
Bluetooth 4.0
automatic pairing

Technical specification:

- » 4th generation diagnostic electronics
- » LED indicators (measurements independent of the application)
- » Bluetooth with on/off option